U.S. EPA, Region 8 Laboratory

5/8/2013

Analytical Method Information

		Reporting	Surrogate	Duplicate	Matrix Spike		Blank Spike / LCS	
Analyte	MDL	Limit	%R	RPD	%R	RPD	%R	RPD
Diesel-Low-1000ml in Water (EPA 8015)								
Preservation: 04-Hydrochloric Acid (HCl) to pH<2, Cool 4°C								
Container: 12_1000mL Amber Glass pH		Amount Required: 1000 mL			Hold Time: 7 days			
<2 w/ HCl								
Diesel range organics	20.0	20.0 ug/L		25	60 - 120	25	60 - 120	25
surr: o-Terphenyl	0.0100		60 - 140	30				